

DRESSING TROLLEY | THR-1025

**DRESSING TROLLEY
WITH BOWL & BUCKET**



DESCRIPTION:

- Approx. Size : 30"x18"x32"
- Frame work made of S.S. tubes
- Two S.S. Shelves over with S.S. railings on all four sides
- The trolley fitted S.S. bucket and S.S. bowl
- The trolley mounted on 100 mm dia castors
- S.S. Frame, S.S. shelves with S.S. Basin & Bucket

| | |
|-------------|----------------------------|
| THR-1025(A) | M.S. FRAME S.S. SHELVES |
| THR-1025(B) | S.S. FRAME S.S. SHELVES |